# TRUE METRIX AIR

Self Monitoring Blood Glucose System

## Quick Reference Guide

### Before You Test

- ✓ Contact your Doctor or Healthcare Professional to determine how often you should test and what your target ranges are for blood glucose results.
- ✓ ALWAYS read the Owner's Booklet before testing for the first time and before using this Quick Reference Guide.

### Meter and Test Strip

Front

Back

## Blood Glucose Testing

**Step 1**  
Wash your hands with warm, soapy water. Dry thoroughly.

**Step 2**  
Gather your materials and prepare your lancing device. Remove test strip from vial and close vial immediately.

**Step 3**  
Insert strip into Test Port with black blocks facing up. The meter turns on.

**Step 4**  
Lance fingertip. Allow blood drop to form. Refer to lancing device Instructions for Use for detailed instructions.

**Step 5**  
Wait until Drop Symbol appears in Display. With strip still in meter, touch Tip of strip to blood drop. Allow blood to be drawn into strip. **Remove test strip from blood drop immediately after meter beeps and dashes appear across Meter Display.**


**Step 6**  
Glucose result will be displayed with alternating date and time.

**Step 7**  
Hold meter in a vertical position with strip pointing down. Press Strip Release Button to discard the used strip from the meter. Testing is complete.

**WARNING:** Upon opening the test strip carton, examine the product for missing, damaged or broken parts. Ensure the test strip vial cap is securely closed. If the product is damaged or the vial cap is not closed, DO NOT use the test strips for testing; product may give inaccurate results. Contact Trividia Health Customer Care at 1-800-803-6025 for replacement and assistance.

### Step 2 Tips


Peel and discard the clear protective sticker from the Meter Display.

Remove one Test Strip from vial. Then immediately press down vial cap firmly to close. See below for additional storage conditions.

### Step 3 Tips

With meter off, insert Test Strip into Test Port with Contact End (black blocks) facing up. The Meter turns on.


**DO NOT place Test Strip in backwards.**

### Step 5 Tips

With strip in Meter, touch Tip of strip to blood drop.

**DO NOT place blood onto Test Strip before inserting Test Strip into Meter Test Port.**

**DO NOT place blood drop on top of Test Strip or any other location except the Sample End (Tip) of the Test Strip.**

### After You Test

- ✓ ALWAYS store the test strips in the original vial with the cap tightly closed.
- ✓ ALWAYS keep the cap of the test strip vial tightly closed.
- ✗ NEVER transfer test strips from one vial to another. Do not expose to heat, moisture, or humidity.
- ✗ NEVER compare your meter to other meters.
- ✗ NEVER immerse meter or test strips in liquids or place any liquids onto meter or into the Test Port.
- ✓ Store TRUE METRIX® AIR System in a dry place at room temperature (40°F-86°F) at 10%-80% relative humidity. DO NOT FREEZE.
- ✓ Dispose of used test strips and lancets in accordance with local regulations.

### Additional Information

- ALL parts of your blood glucose monitoring system could carry blood-borne pathogens after use, even after cleaning and disinfecting.<sup>1,2</sup>
- Cleaning and disinfecting the lancing device and the meter destroys most, but not necessarily all, blood-borne pathogens.
- Wash your hands thoroughly with soap and warm water before and after handling the meter, lancing device, lancets, or test strips as contact with blood presents an infection risk.
- If the meter is being operated by a second person who provides testing assistance, the meter and lancing device should be disinfected prior to use by the second person.
- It is important to keep the meter and the lancing device clean and disinfected. For instructions on how to clean and disinfect the meter and lancing device, see Meter Cleaning and Disinfecting and Lancing Device Cleaning and Disinfecting in the Owner's Booklet.

Customer Consultation  
English or Spanish, Monday - Friday,  
8AM - 8PM EST, 1-800-803-6025  
[www.trividiahealth.com](http://www.trividiahealth.com)

Manufactured by:

**TRIVIDIA**  
HEALTH™

FORT LAUDERDALE, FL 33309 U.S.A.

1. FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication. Available at <https://wayback.archive-it.org/7993/2017011013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>. 2. Infection Prevention during Blood Glucose Monitoring and Insulin Administration. Available at <http://www.cdc.gov/injectionsafety/blood-glucose-monitoring.html>.

© 2023 Trividia Health, Inc. TRUE METRIX and Trividia Health logo are trademarks of Trividia Health, Inc. REA4TVH08 Rev. 43

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